

Cyres Cytology – Upgrade Features

The latest version of Cyres Cytology – version 4.6 – is now available. This includes a range of new features and we recommend all Labs to upgrade to this version as soon as possible – your system will check for upgrades automatically and try and download them from the web. If your firewall prevents this please ask us for an upgrade CD. New features include:

1. A wide range of new turnaround reports to help Labs achieve the new 14 day turnaround time. These help identify where bottlenecks or delays may be occurring by measuring turnaround times between primary screeners / QC screeners / checkers / consultants etc. There are also all the reports showing date sample taken to date sample received and on to date sample authorised.
2. APV and PPV reports for checkers.
3. New reasons for inadequacy for LBC samples – in line with comments via the Cyres Blog these reasons have been updated within a new drop down option. The old reasons have been kept but moved to a new ‘PAP’ tab.
4. The system now has a new access level – Failsafe Officer – which enables the Failsafe module but nothing else.
5. There are also enhancements to failsafe and system administrator options as well as improved support for scrolling mice.
6. All of the above are of course in addition to all the new features contained within previous upgrades including support for SQL Server (this dramatically improves performance so anyone who has yet to take up this option is urged to do so), enhanced staff reports and so on.

We look forward to hearing back from you – should you need any further information or help with upgrading please do not hesitate to contact us on the numbers below.

